

OMEGA TYPE ES-316 Experimental Set Up has been designed specifically for determine wave length of Laser Light by diffraction grating method. The set-up consists of Circular Table, Diode Laser, Laser detector, Diffraction Grating, Nanoammeter, Reading Lens and Spirit Level.

The set up is complete in all respect and requires no other apparatus. Practical experience on this set up carries great educative value for Science and Engineering Students.

OBJECT

To determine wave length of Laser Light by diffraction grating method.

FEATURES

The complete Experimental Set-up consists of the following items.

01. Diode Laser with Power Supply.

- Maximum output : 2 mW
- Wave length : 670 nm visible red
- Power supply : Included with ON/OFF switch working on 230 mains.

02. Circular Table : Spectrometer scale 6" dia circle with vernier but without Collimator & Telescope. It has two holders one for laser & other for Laser detector.

03. Laser Detector : Composition silicon Laser detector mounted in Aluminium case.

04. Diffraction Grating : Hilger & Watts Type, 15000 lines per inch/6000 lines per cm.

05. Nanoammeter : OMEGA TYPE DNM-021

06. Reading Lens : 50 mm diameter with handle

07. Spirit Level : 60 mm length

* Weight : 9.2 Kg. (Approx.)

* Strongly supported by detailed Operating Instructions, giving details of Object, Theory, Design procedures, Report Suggestions and Book References.

OMEGA

We are committed to the continuous development of our products, and therefore reserve the right to amend specifications without prior notice.

OMEGA ELECTRONICS

India Sales

B-28, Fateh Singh Scheme , Opp. Rajputana
Palace Sheraton, Jaipur-302006, INDIA
Ph: +141-2375647 / 2379223
Fax: +91-141-2204481

Email: omega@sancharnet.in;
info@omegaelectronics.net
Web: www.omegaelectronics.net

Mfg./ Export/ Training

28 E & F, Malviya Industrial Area,
Jaipur-302017, INDIA
Ph: +91-141-2751136 / Fax: +91-141-2751559